

# Fortera<sup>®</sup>

*Fortifying Osteoconductive Biocement*



- Injectable, self-setting calcium phosphate cement with superior compressive strength


## Performance Characteristics

Composition	Fully synthetic
Compressive Strength	45 MPa
Durable Strength	Unlike cements derived from calcium sulfate, Fortera will not wash away or rapidly degrade
Injectable	16 gauge
Osteoconductive	Material crystallizes into hydroxyapatite upon setting, the mineral phase of native human bone
Thermally Stable	Does not require refrigeration prior to use
Isothermic setting	Safe for adjacent tissue contact
Resistant to lavage	Solidifies in a wet field

### Products

Description	Kit Size
Fortera	3cc
Fortera	5cc
Fortera	10cc

### Superior Compressive Strength



Fortera exhibits superior compressive strength when compared to industry leading calcium-based cements<sup>†</sup>

#### Vivorté, Inc.

1044 East Chestnut Street  
Louisville, Kentucky 40204